



AltaFlex® is the next generation of CPAP tubing featuring FlexFilm™ Technology



ELIMINATES
Stiffness and bend memory.



MOVES EFFORTLESSLY
AltaFlex® tubing moves effortlessly during sleep to minimize mask leaks.



THIS MEANS...
Better Therapy and Better Sleep!



Certification EN ISO 13485:2016 Medical Device



FDA Registered



CE Mark



UDI Marking Compliance 21 CFR 801.45

CONTACT US



+(1) 949-581-9530

Get in **Touch**
+(1) 949-581-9530

For more info, visit our **www.AltaFlex.US**
Email: sales@smoothborplastics.com



Smooth-Bor[®] has manufactured tubing in the USA for 40+ years. Our business was built on innovation and design partnerships with the largest CPAP manufacturers in the world.



Advanced Technology

- 100% Made in the USA
- 15mm & 19mm Internal Diameters
- Standard ISO 22mm Cuff Connector
- Length: 72 inches
- Full Biocompatibility Testing to 10993 & 18562
- UDI Marked: Compliance to 21 CFR 801.45
- LDPE & TPE Prime Materials
- Ergonomic Cuff Design
- Fits Most Major Brands of CPAP Machines

Get in **Touch**
+(1) 949-581-9530

For more info, visit our **www.AltaFlex.US**
Email: sales@smoothborplastics.com